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THE HEAT EMERGENCY KIT HELPS SAVE LIVES

CONTENT LIST

- 1 x E.COOLINE Repower cooling blanket
- 0.5 L of activation water for the cooling blanket
- 2 Paper towels
- 1 Forehead thermometer card
- 1 Pulse Oximeter (not available in basic version)
- 1 Drinking cup
- 2 x 100 ml Drinking water
- 2 Sachets of electrolyte plus powder, 6 g each
(to dissolve in the cup with drinking water from bag)
- 2 x Medical face mask
- 2 x Resuscitator Face Mask / CPR face shield
- 2 pairs of Disposable gloves Showa
- 2 x Spit bags
- 2 x Disinfectant wipes Sterillium Tissue

ATTENTION: Please check the expiration date of the individual products regularly and replace them in good time. Please also remember to check that the batteries are working properly.

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Patenterte Technologie | Patented Technology

HITZENOTFALL-SERVICE HEAT EMERGENCY SERVICE

pervormance international GmbH

Mühlsteige 13 · 89075 Ulm

Tel. +49 731 14 071-0 · Fax +49 731 14 071-14

info@e-cooline.de · www.e-cooline.de

HITZENOTFALL? HEAT STRESS?

HITZENOTFALL-KOFFER HEAT EMERGENCY KIT

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Made in Germany
Zum Patent angemeldet
Patent pending

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For cooling & aid, please read instruction

pervormance international GmbH
Mühlsteige 13, 89075 Ulm, Germany
www.e-cooline.de

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Heat doesn't only affect vulnerable groups such as older people, pregnant women and children. But also other people are suffering from heat problems and are at risk due to high temperatures. According to the German Weather Service, heat waves are becoming more and more common in summer.

Even in otherwise rather cool countries such as Canada, temperatures of 40-50 °C have already been registered. That means sweating a lot. But cooling down the body through sweat puts a strain on the body. Over 30,000 deaths in recent years were due to heat. Heat also caused 108,000 hospital admissions and emergencies in 2020: 177% more than in the past.

The number of deaths has increased more than eightfold. Cities, municipalities and companies are therefore already developing heat action plans. Protect yourself too!

YOU CAN RECOGNIZE HEAT COLLAPSE BY THE FOLLOWING SYMPTOMS:

- High body temperature up to 38 – 40 °C
- Cool, pale, moist skin
- Heavy sweating
- Fast but weak pulse
- Lightheadedness, headaches and nausea
- Weakness, dizziness

YOU CAN RECOGNIZE A HEAT STROKE BY THE FOLLOWING SYMPTOMS:

- High body temperature up to 40 °C or more
- Hot, dry and red skin
- No sweating!
- Fast, strong pulse
- Dizziness, headache and vomiting
- Muscle cramps
- Drowsiness and possibly unconsciousness

If a heat stroke is not noticed and counteracted in time, it can result in:

- a drop in blood pressure,
- circulatory collapse and associated shock
- as well as multiple organ failure.

Due to the collapse of the entire cardiovascular system, the organs are no longer adequately supplied with blood and oxygen. This can lead to, for example, shortness of breath, fainting or seizures. Without cooling and medical treatment, heat stroke can, in the worst case scenario, lead to death.

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FIRST AID MEASURES:

- Move the person into a shady area.
- Take off constricting or excessively warm clothing.
- Check their consciousness and breathing until the emergency service arrives.
- If the affected person is fully conscious and sitting upright, dissolve electrolytes in the drinking cup with drinking water and offer it to them.
- Temperature check with forehead thermometer card.
- Clip the pulse oximeter onto a finger, switch it on and measure vital signs
- If the person vomits, give them the spit bag and use the enclosed hygiene products, such as disinfectant wipes, mask and gloves.

FOR EXPERIENCED FIRST AIDERS:

- If the person in need of help is not breathing, place the enclosed respiratory mask over their mouth with the valve facing upwards. Secure the two rubber bands on the sides of the ventilation sheet behind their ears. Ventilate the person through the valve and keep checking their breathing.
- If the person is unconscious, perform chest compressions. To help you find the perfect rhythm for the chest compressions of 100 to 120 beats per minute, you can use the rhythm of the songs „Stayin' Alive“ by the BeeGees or „Highway to Hell“ by AC/DC as a guide.

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STEP 1: DIAL EMERGENCY NUMBER

STEP 2: COOLING IN HEAT EMERGENCIES – HOW IT WORKS

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IMPORTANT NOTE: The E.COOLINE Repower cooling blanket is only intended for service first aid in the event of an acute heat emergency and does not replace medical examination or treatment. Please contact an emergency doctor/emergency service. The contents of the heat emergency kit are only intended to be used for first aid until the doctor/emergency services arrives.