Dräger DrugCheck 3000

Instructions for use

de, enUS, fr, es, ptBR

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1 Safety-related information

1.1 Basic safety instructions

– Comply with all local and national rules and regulations associated with this product.

– Notify Dräger in the event of any component fault or failure.

– Do not use a faulty or incomplete product. Do not modify the product.

– Only trained and competent users are permitted to use this product.

– Do not dispose of the instructions for use. Ensure that they are stored and used appropriately by the user.

– Strictly follow the instructions for use. The user must fully understand and strictly observe the instructions. Use the product only for the purposes specified in the intended use section of this document.

1.2 Precautions and warnings

1.1 To ensure proper hygiene, use gloves when handling the Dräger DrugCheck 3000. Do not touch the white end of the sample collector before or after sampling. Observe the common rules of hygiene.

– Follow proper handling and disposal procedures.

– Open the foil pouch only immediately before use to prevent contamination of the sample collector.

– Do not use the Dräger DrugCheck 3000 if the foil pouch is damaged in any way.

– Do not use the Dräger DrugCheck 3000 beyond the expiration date indicated on the foil pouch. The expiry date is in the format YYYY-MM. Example: 2024-01 means that the Dräger DrugCheck 3000 is not permitted to be used after the end of January 2024.

2 Description

2.1 Product overview

2.1.1 Dräger DrugCheck 3000 Test Kit (Fig. A)

1 Sample collector with color indicator
2 Test cassette with indicator ring
3 Funnel-shaped test opening
4 Security seal
5 Result window with test strips

6 Control lines area (C)
7 Test lines area (MET, AMP, THC, BZO, OPI, COC)

Each packaged unit contains:

– 20 Dräger DrugCheck 3000 test kits, each individually packaged in a sealed foil pouch
– 1 x Instructions for Use

The Dräger DrugCheck 3000 test kit consists of 2 components: the test cassette and the sample collector.

The removable sample collector has a discolorable red sample adequacy indicator.

The test cassette contains:

– One vial containing buffer fluid
– A red indicator ring for handling control
– 2 test strips with 2 control lines (C)
– up to 6 test lines for the different drugs or drug classes, depending on the product variant

Additional recommended material:

– Disposable gloves (e.g., latex or nitrile gloves)
– Timer

2.2 Intended use

1.1 Not intended for medical use.

The Dräger DrugCheck 3000 is intended for forensic applications as well as for use by laypersons and in the workplace. Observe local and national regulations for any of the described intended uses.

The Dräger DrugCheck 3000 is designed for simultaneous, qualitative detection of up to 6 drugs or drug types (depending on the product variant) in human saliva samples. These drugs or drug types are cocaine, opiates, amphetamines, methamphetamines, THC/cannabis, and benzodiazepines. Each drug or drug type is represented by a separate test line on one of the test strips in the result window on the test cassette.

The Dräger DrugCheck 3000 packaging gives all key information on the drugs or drug types with their abbreviations (see also 2.2.1), detection limits (in ng/mL), shelf life, and storage conditions.

The detection limit (cut-off) is the deciding factor as regards the presence of the drugs or drug types in the sample. The test result is correspondingly positive or negative above or below the limit value. To obtain a confirmed test result, a more specific method must be used. The preferred methods are liquid chromatography or mass spectrometry coupled with mass spectrometry (LC-MS or GC-MS). All the results from the Dräger DrugCheck 3000 require professional assessment with reference to further clinical assessment of the test subjects. This applies in particular to a positive („non-negative“) test result.

2.2.1 Abbreviations

Droge	Target analyte
AMP	Amphetamine
BZO	Benzodiazepine
COC	Cocaine
MET	Methamphetamine
OPI	Opiates
THC	Tetrahydrocannabinol
	Δ9-Tetrahydrocannabinol

2.3 Explanation of symbols

Symbol	Explanation
1	Follow the Instructions for Use.
2	Manufacturer

Symbol	Explanation
2	Intended for single use only
LOT	Batch
REF	Order number
5	Temperature range for storage conditions
6	Expiration date
7	Contents sufficient for <= 6 tests
8	Measurement mode: fast (THC)
9	Measurement mode: sensitive (THC)

2.4 Test principle
The Dräger DrugCheck 3000 is based on the immunoassay principle of competitive inhibition. Drugs contained in the sample compete with drugs on the test membrane for the cohesion of nanoparticles coated with antibodies.

Sampling:
The Dräger DrugCheck 3000 is designed for use with oral fluid samples, which are taken with the integrated sample collector.

The sample does not require special treatment. The sample is taken by taking a drop of the sample collector by absorption. The sample collector is then inserted into the funnel-shaped test opening in the test cassette to perform the test.

Test:
The oral fluid sampler is then rinsed in the buffer solution by shaking. After waiting for the pre-incubation time to elapse, the test can be started by breaking the security seal and pressing the sample collector downwards with the test cassette.

Gold particles coated with antibodies react with drug molecules from the fluid sample and, during the course of the test, with drug conjugates on the test membrane. If the sample is drug free, the antibodies are free to react with the drug conjugates. When this happens, you will see a red line on the test strip. If the drug is available in a sufficient concentration, the gold particles are impeded from adhering to the drug conjugate. If the sample is preliminarily positive („not negative“), no test line will appear in the area for the corresponding drug in the test strip result window.

Quality control:
An additional function check independent of the fluid sample is integrated into every test membrane (Fig. B/6). The appearance of a red line in area C in the viewing window indicates that sufficient saliva was collected and the test membrane fully wetted and that the antibody-antigen-reaction occurred as expected.

en US
Instructions for use

1 Safety-related information

1.1 Basic safety instructions

– Results may be incorrect if storage times are not properly observed.

3.1 Preparing the test

Prerequisites

– Adhere to the specified storage times.

– Ensure that the test subject has not eaten or drunk anything for at least 10 minutes before sampling, e.g., chewing gum, tobacco, coffee, or other substances.

– Ensure that the ambient temperature is between +5 °C to +40 °C in a non-condensing environment.

– Open the foil pouch.

– Only remove the sample collector immediately before using the test kit.

3.2 Carrying out the test

Collecting the sample

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© Dräger Safety AG & Co. KGaA