

Chem-Tape® Spezialklebeband | Art.- Nr. CT54

Artikelbeschreibung: 55-Meter-Rolle, 5 cm breit, ausgezeichnete Chemikalienbarriere, zum Abkleben von Handschuhen und Overalls, geprüft auf den Durchbruch diverser aggressiver Chemikalien. Die notwendige Ergänzung im Einsatz von Chemikalienschutzanzügen und Handschuhen.


Farbe: gelb

VE: 24 Rollen

Menge pro Palette: 720 Stück/ 30 Kartons

Kartonmaße/ cm: ca. 35x35x30,5

Kartongewicht/ kg: ca. 23



Chem-Tape®			
Chemical	Breakthrough Time	Permeation Rate	SDL (ppm)
	Normalized¹ (min)	(ug/cm²/min)	
Acetone	>480	ND	0.020
Acetonitrile	>480	ND	0.020
Acrylonitrile	>480	ND	0.190
Ammonia Gas	>480	ND	0.500
Ammonia Liquid	>480	ND	0.500
Benzene	123	23	0.100
Carbon Disulfide	>480	ND	0.100
Chlorine Gas	>480	ND	0.050
Chlorine Liquid	>480	ND	0.050
Dichloromethane	>480	ND	0.020
Diethylamine	167	7.22	0.020
Dimethylformamide N,N-	>480	ND	0.020
Ethyl Acetate	>480	ND	0.020
n-Hexane	>480	ND	0.020
Hydrofluoric Acid (48-51%)	>480	ND	0.100
Hydrogen Fluoride Anhydrous	>480	ND	0.200
Methanol	>480	ND	0.020
Mustard (HD) Chemical Agent	>480 ²	NA	NA
Nerve (VX) Chemical Agent	>480 ²	NA	NA
Nitric Acid (70%)	>480	ND	0.025
Nitrobenzene	>480	ND	0.020
Sarin (GB) Chemical Agent	>480 ²	NA	NA
Sodium Hydroxide (50%)	>480	ND	0.050
Soman (GD) Chemical Agent	>480 ²	NA	NA
Sulfuric Acid	>480	ND	0.050
Tetrachloroethylene	>480	ND	0.020
Tetrahydrofuran	>480	ND	0.020
Toluene	>480	ND	0.020

ND = None Detected NA = Not Applicable

¹ Breakthrough Times normalised to a permeation rate of 0.1µg/cm²/min in accordance with ASTM F739

² Testing conducted per MIL-STD 282 and US Army Test Operations Procedure (TOP) 8-2-501

Sources for all test data are independent laboratories. All tests were performed under laboratory conditions and not actual use conditions.

Wayne Harwood
Technical Support Manager

